

Wall Grab Bars and Car Handle Assist for Elderly & Handicap Users

Introduction:

Welcome to our wall grab bars and car handle assist system designed with the utmost care and attention to detail for elderly and handicap individuals. This user manual aims to guide you through the installation process, safe usage, and maintenance of your new grab bars, ensuring optimal safety and comfort.

Safety Precautions:

- *Always consult a healthcare professional or occupational therapist before installing grab bars in your home or vehicle.*
- *Ensure the installation surface (wall, door frame, or car door) is strong enough to support the grab bars' weight and intended usage.*
- *Avoid using the grab bars for purposes other than providing support or assistance.*

Installation Guide:

1. **Determine Location:** *Choose a strategic location where the grab bars will be most useful, such as near doorways, stairs, or in the car.*
2. **Prepare Surface:** *Ensure the installation surface is clean, dry, and free of obstructions.*
3. **Assemble Components:** *Verify that all components (grab bars, screws, anchors, if applicable) are included and in good condition.*
4. **Mark & Drill:** *Using a pencil, mark the drill points on the installation surface according to the pre-drilled holes in the grab bars. Drill pilot holes if necessary.*

5. **Install:** Position the grab bars and secure them to the wall or door frame using the provided screws. Tighten the screws evenly to ensure a firm hold.
6. **Test & Adjust:** Once installed, test the grab bars for stability and adjust as needed.

Usage Instructions:

- Grip the grab bars firmly, using the grip grooves for added stability.
- Do not hang or pull on the grab bars excessively; they are designed for support during transitions, not for weight-bearing activities.
- In vehicles, use the car handle assist carefully and follow the manufacturer's guidelines for your specific vehicle model.

Maintenance & Cleaning:

- Regularly inspect the grab bars for signs of wear, damage, or loosening screws.
- Clean the grab bars with a damp cloth and mild soap or disinfectant, avoiding abrasive chemicals that may damage the surface.
- Do not use excessive force or sharp objects when cleaning to prevent scratching or damaging the non-slip surface.

Warranty & Support:

- Please refer to the warranty card included with your purchase for specific warranty information and contact details.

- *If you encounter any issues or have questions during installation or usage, feel free to contact our customer support team for assistance.*

Thank you for choosing our grab bars and car handle assist system. We hope they contribute to a safer and more independent lifestyle for you or your loved ones.