

Shower Seat Use and Care Manual

This shower seat is designed for durability and ease of use in wet environments, featuring a corrosion-resistant 304-grade stainless steel frame and a waterproof phenolic top. It is rated to support up to 400 pounds when properly installed, meeting ADA guidelines for accessible design, suitable for both home and commercial use.

Before using the shower seat, ensure it is properly installed and locked securely into place. To deploy the seat, gently unfold it until it clicks into a stable, horizontal position; do not use the seat unless it is fully locked to avoid slipping or accidents.

When using the seat, sit centrally on the phenolic top to distribute weight evenly. Do not exceed the maximum weight capacity of 400 pounds, as this may cause damage to the seat or lead to safety hazards. Avoid sudden movements or standing on the seat, as it is designed for seating only.

For users with mobility challenges or seniors, use the seat slowly and carefully. Ensure there is no water accumulation on the seat surface before sitting to prevent slipping; the phenolic top is waterproof, but excess water should be wiped away for added safety.

After use, fold the seat smoothly against the wall to save space and keep the shower area accessible. Ensure the seat is fully folded and secured to avoid it accidentally unfolding when not in use.

To maintain the shower seat's performance and appearance, clean it regularly. The waterproof phenolic top can be cleaned quickly with a damp cloth and mild soap; avoid using harsh chemicals, abrasives, or scouring pads, as these may scratch or damage the surface.

Wipe the 304-grade stainless steel frame with a dry cloth after each use to prevent water spots and corrosion. For stubborn stains, use a non-abrasive cleaner suitable for stainless steel, then rinse and dry thoroughly.

Check the seat regularly for loose parts, including mounting hardware, frame connections, and the locking mechanism. If any parts are loose or damaged, stop using the seat immediately and tighten or replace the parts according to the manufacturer's recommendations.

Avoid exposing the seat to extreme temperatures or harsh chemicals, as this may compromise the durability of the materials. Do not use the seat if it shows signs of damage, such as cracks, rust, or loose locking mechanisms.

Proper care and use will extend the service life of the shower seat. If you have any questions about use or maintenance, refer to the manufacturer's detailed instructions or seek professional assistance.