



up to **\$40,000**

Excavation/

Replacement



Reclaim Your Home and Yard!

Save your septic system with the patented and proven power of Aerobic Action.

Eliminate foul septic odor and pooling liquid in septic drain field—prevent sewage backflow and back up.

up to

\$7,000

Chemicals/

Additives

Fix sluggish drains, gurgling pipes and toilets that are slow to flush.

Save Thousands

of Dollars!

Costs Under

\$1,200

SAVER

Your septic tank discharge will be 90% cleaner than before, a benefit to your local groundwater.

University Studies Agree

"... adding an aerobic treatment device restores the original functioning of the soil absorption system (drainfield) in a matter of months..."

— Purdue University

"... after nine months of operation, the aerobic systems were discharging effluent 95% cleaner than before ... clean enough to surface discharge ..." — Baylor University

Photograph the code with your smart phone! www.septicsystemsaver.net



12" x 14" x 10" • Weather proof • Plugs into 115V outlet

Five Year Cost Comparison

SepticSystemSaver.net



Drill 5/8" hole for airline to exit septic tank



Dig trench and bury airline 2"- 3" underground



Model: S3-100 UPC: 793573194879 Septic System Saver®

- treatment capacity 4 people max.
- 6' power cord
- 20' airline
- 100 Watts 0.75 amp
- \$5 \$7 per month to operate
- 1 meter of Bio-Brush[™]
- quiet operation 39 dB
- 12 month warranty

GOOD Restoration



Model: S3-101 UPC: 793573194862 Septic System Saver®-Classic

- treatment capacity 6 people max.
- 6' power cord
- 20' airline
- 115 Watts 1.0 amp
- \$6 \$8 per month to operate
- 2 meters of Bio-Brush[™]
- quiet operation 39 dB
- 12 month warranty

Better Restoration





Plug into standard outlet and watch it work



(800) 969-8795

in Four Steps!

Model: S3-102 UPC: 793573062994 Septic System Saver®-PRO

- treatment capacity 8 people max.
- 6' power cord
- 20' airline
- 150 Watts 1.3 amp
- \$7 \$10 per month to operate
- 4 meters of Bio-Brush[™]
- guiet operation 39 dB
- 24 month warranty

BEST Restoration



Note: Selecting a higher capacity model than required will speed up the restoration process! Selecting too small of a capacity will hinder the restoration process.